

# Experiment #1 Laser Beam Parameters


EP20BC015

Sameeth Kishore (Guide)

Aim To determine the various parameters of a Laser Beam, specifically Divergence ( $\theta_0$ ), Beam Radius ( $w_0$ ), and Location of Beam Radius ( $z$ )

Apparatus Laser Beam Generator, meter scale, detector, knife blade, Photo detector, lens, Power Meter, etc

Diagram



Theory The technique that will be used to determine Power distribution, spot size, divergence, etc is as follows:  
We shine the laser in the direction of a knife edge that is slowly inserted into the path of the laser. Further, the laser is assumed to be oscillating in TEM<sub>00</sub> mode so that the spatial distribution of the beam is Gaussian. Hence, the irradiance distribution  $I(x, y)$  is

$$I(x, y) = \frac{2P_0}{\pi w_0^2} \exp\left(-\frac{2(x^2 + y^2)}{w_0^2}\right)$$

Where  $P_0$  Total Power of the beam  
 $w_0$  (Spot size of beam) / 2

Power transmitted past a knife edge blocking of all points  $x < a$

$$\begin{aligned} P &= \int_{y=-\infty}^{\infty} \int_{x=a}^{\infty} I(x, y) dx dy = \frac{2P_0}{\pi w_0^2} \int_{y=-\infty}^{\infty} \int_{x=a}^{\infty} \exp\left(-\frac{2(x^2 + y^2)}{w_0^2}\right) dx dy \\ &= \frac{P_0}{2} \operatorname{erfc}\left(\frac{a\sqrt{2}}{w_0}\right) \end{aligned}$$

Now, 
$$\frac{I(x,y)}{P_0} = \frac{1}{2\pi(w')^2} \exp\left(-\frac{(x^2+y^2)}{2w'^2}\right)$$
  $w'$  - standard deviation

Normalized Gaussian Form

and 
$$\frac{P}{P_0} = \frac{1}{2} \operatorname{erfc}\left(\frac{a}{\sqrt{2} w'}\right)$$

It is easy to show that the points with 25% and 75% of relative powers are located at distances equal to the probable error  $e_p$  ( $e_p \approx 0.6745 w'$ ) on either side of the Gaussian Distribution's maximum. Hence,  $w'$  can be obtained from the relative power vs. knife-edge position curve and beam spot size ( $2w_0 = 4w'$ ) can be calculated.

To determine the Divergence of the laser beam, we once again assume that the laser is oscillating in  $TEM_{00}$  mode, and its output is a Gaussian beam as before. It can be shown that:

$$w(z) = w_0 \left( 1 + \left( \frac{\lambda z}{\pi w_0^2} \right)^2 \right)^{1/2}$$

where  $w(z)$ : Radius of the beam

$z$ : axis along the direction of propagation

$\lambda$ : wavelength of laser radiation

$w_0$ : (spot size of beam) / 2.

The gradient  $\frac{dw(z)}{dz} = \theta(z)$  of the beam radius at  $z$  is:

$$\theta(z) = \left( \frac{\lambda}{\pi w_0} \right)^2 \frac{z}{w_0} \left( 1 + \left( \frac{\lambda z}{\pi w_0^2} \right)^2 \right)^{1/2}$$

when  $z \rightarrow \infty$ , ( $\theta(z) \rightarrow \theta_0$ ) where  $\theta_0 = \frac{\lambda}{\pi w_0} = 1/2$  Divergence

$w^2(z) = w_0^2 + \theta_0^2 z^2$ . Here,  $w_0$ ,  $\theta_0$ , and  $z$  are all unknown

we choose 3  $z$  values:  $z$ ,  $z+D$ ,  $z+2D$  so that we get:

$$\left. \begin{aligned} w_1^2 &= w_0^2 + \theta_0^2 z^2 \\ w_2^2 &= w_0^2 + \theta_0^2 (z+D)^2 \\ w_3^2 &= w_0^2 + \theta_0^2 (z+2D)^2 \end{aligned} \right\} \begin{aligned} &\text{Form a set of equations which} \\ &\text{we can solve to find} \\ &w_0, \theta_0, \text{ and } z. \end{aligned}$$

From these equations, we obtain:

$$2\theta_0^2 D^2 = \omega_3^2 - 2\omega_2^2 + \omega_1^2 \Rightarrow \theta_0 = \frac{1}{\sqrt{2}D} (\omega_3^2 - 2\omega_2^2 + \omega_1^2)^{1/2}$$

from this  $\theta_0$ ,  $\omega_0 = \frac{\lambda}{\pi\theta_0}$  and  $z = \frac{(\omega_1^2 - \omega_0^2)^{1/2}}{\theta_0}$

In this way, Divergence, Beam Spot Size, and Location of Beam Radius, can be calculated.

### Observations:

(i) Set of measurements with  $z = 300\text{mm}$ : (To determine  $\omega_1$ ):

| S. No. | Position of knife (mm) | Power Meter Reading (mW) | Normalized |
|--------|------------------------|--------------------------|------------|
| 1.     | 3.4                    | 4.97                     | 1.0        |
| 2.     | 3.3                    | 4.90                     | 0.986      |
| 3.     | 3.2                    | 4.75                     | 0.956      |
| 4.     | 3.1                    | 4.34                     | 0.873      |
| 5.     | 3.0                    | 3.78                     | 0.761      |
| 6.     | 2.9                    | 2.88                     | 0.579      |
| 7.     | 2.8                    | 1.90                     | 0.382      |
| 8.     | 2.7                    | 1.09                     | 0.219      |
| 9.     | 2.6                    | 0.057                    | 0.011      |
| 10.    | 2.5                    | 0.027                    | 0.005      |
| 11.    | 2.4                    | 0.010                    | 0.002      |
| 12.    | 2.3                    | 0.003                    | 0.0006     |
| 13.    | 2.2                    | 0.001                    | 0.0002     |
| 14.    | 2.1                    | 0                        | 0          |

(ii) Set of readings with  $z = z + D = 400\text{mm}$ : (To determine  $w_2$ )

| S. No. | Position of Knife (mm) | Power-Meter Reading (mw) | Normalized. |
|--------|------------------------|--------------------------|-------------|
| 1.     | 6.2                    | 4.94                     | 1.0         |
| 2.     | 6.1                    | 4.86                     | 0.984       |
| 3.     | 6.0                    | 4.58                     | 0.927       |
| 4.     | 5.9                    | 4.10                     | 0.830       |
| 5.     | 5.8                    | 3.34                     | 0.676       |
| 6.     | 5.7                    | 2.51                     | 0.508       |
| 7.     | 5.6                    | 1.69                     | 0.342       |
| 8.     | 5.5                    | 0.97                     | 0.196       |
| 9.     | 5.4                    | 0.46                     | 0.093       |
| 10.    | 5.3                    | 0.17                     | 0.034       |
| 11.    | 5.2                    | 0.066                    | 0.013       |
| 12.    | 5.1                    | 0.030                    | 0.0061      |
| 13.    | 5.0                    | 0.015                    | 0.003       |
| 14.    | 4.9                    | 0.010                    | 0.002       |
| 15.    | 4.8                    | 0.0                      | 0           |



1058  $\frac{1}{2}$

(iii) Set of readings with  $z = z + 2D = 500 \text{ mm}$  (To determine  $w_3$ )

| S.No. | Position of knife (mm) | Power-Meter Reading (mw) | Normalized |
|-------|------------------------|--------------------------|------------|
| 1.    | 11.4                   | 4.96                     | 1.0        |
| 2.    | 11.3                   | 4.89                     | 0.986      |
| 3.    | 11.2                   | 4.70                     | 0.948      |
| 4.    | 11.1                   | 4.34                     | 0.875      |
| 5.    | 11.0                   | 3.75                     | 0.756      |
| 6.    | 10.9                   | 2.91                     | 0.588      |
| 7.    | 10.8                   | 2.19                     | 0.441      |
| 8.    | 10.7                   | 1.42                     | 0.286      |
| 9.    | 10.6                   | 0.93                     | 0.187      |
| 10.   | 10.5                   | 0.44                     | 0.090      |
| 11.   | 10.4                   | 0.14                     | 0.028      |
| 12.   | 10.3                   | 0.05                     | 0.010      |
| 13.   | 10.2                   | 0.02                     | 0.004      |
| 14.   | 10.1                   | 0.01                     | 0.002      |
| 15.   | 10.0                   | 0.005                    | 0.001      |
| 16.   | 9.9                    | 0.0                      | 0.0        |

### Calculations:

(i) For  $z = 300 \text{ mm}$ : 75% power = 0.405 mm  
25% power = 0.681 mm

$$\therefore 2 \times \text{error} = |75\% \text{ power} - 25\% \text{ power}| = |0.405 - 0.681|$$
$$\Rightarrow \text{error} \approx \frac{0.276}{2} \approx 0.138$$

$$\therefore w_1 = \frac{0.138}{0.6745} \approx 0.205 \text{ mm}$$

(300mm)

(ii) For  $z = 400 \text{ mm}$ : 75% power = 0.350 mm  
25% power = 0.660 mm

$$\therefore 2 \times \text{error} = |75\% \text{ power} - 25\% \text{ power}| = |0.350 - 0.660|$$
$$\Rightarrow \text{error} \approx \frac{0.31}{2} \approx 0.155$$

$$\therefore w_2 = \frac{0.155}{0.6745} \approx 0.230 \text{ mm}$$

(400mm)

(iii) For  $z = 500 \text{ mm}$ : 75% power = 0.395 mm  
25% power = 0.744 mm

$$\therefore 2 \times \text{error} = |75\% \text{ power} - 25\% \text{ power}| = |0.395 - 0.744|$$
$$\Rightarrow \text{error} \approx \frac{0.349}{2} \approx 0.1745$$

$$\therefore w_3 = \frac{0.1745}{0.6745} \approx 0.259 \text{ mm}$$

(500mm)

Now,  $\theta_0$  (Divergence of ~~light~~ laser beam) is given by:

$$\theta_0 = \frac{1}{\sqrt{2D}} (w_3^2 - 2w_2^2 + w_1^2)^{1/2} \quad \text{where } D = 100 \text{ mm}$$

$$= \frac{1}{\sqrt{2(100)}} ((0.259)^2 - 2(0.23)^2 + (0.205)^2)^{1/2}$$

$$\therefore \theta_0 = \frac{1}{\sqrt{200}} (0.0671 - 0.1058 + 0.042)^{1/2}$$

$$\approx \frac{1}{\sqrt{200}} \sqrt{0.003325} \approx \sqrt{1.6625 \times 10^{-5}} \approx 0.0041.$$

$$\therefore \theta_0 \approx 4.1 \times 10^{-3}. \text{ Divergence of Laser Beam.}$$

Now,  $w_0 = \frac{\lambda}{\pi \theta_0}$  where  $\lambda$ : wavelength of the laser used  
(He-Ne lasers are of  $\lambda = 632.8 \text{ nm}$ ).

$$\therefore w_0 \approx \frac{632.8 \text{ nm}}{(3.1416)(4.1)(10^{-3})} \approx 4.913 \times 10^4 \text{ nm} \approx 49.13 \mu\text{m}.$$

and finally,  $z = \frac{(w_1^2 - w_0^2)^{1/2}}{\theta} = \frac{((0.205)^2 - (0.049)^2)^{1/2}}{4.1 \times 10^{-3}}$

$$z \approx \frac{0.199}{4.1 \times 10^{-3}} \approx 48.537 \text{ mm}.$$

### Error Analysis:

We measure  $w_1, w_2$ , and  $w_3$  with the help of position of knife edge (with L.C = 0.02 mm) and Power-Meter readings (with L.C = 0.001 mW).

and  $(0.6745) w(2) = P_{25\%} - P_{75\%}$   $P_{25\%}$ : position with 25% Power.

$$\Rightarrow (0.6745)(2) dw = dP_{25\%} + dP_{75\%}$$

Now,  $\text{Power}_{25\%} = \text{Power}_{100\%} \left(\frac{1}{2}\right) \text{erfc}\left(\frac{P_{25\%}}{\sqrt{2}w}\right)$

$$\therefore d\text{Power}_{25\%} = \text{Power}_{100\%} \left(\frac{1}{2}\right) \left(-\frac{2}{\sqrt{\pi}} + \right)$$

Error in measuring  $w$  is given like so:

$$(0.6745) w(2) = x_{25\%} - x_{75\%}$$

$$\Rightarrow 2(0.6745) dw = 2 dx \Rightarrow dw = \frac{1}{0.6745} dx$$

$$\therefore \frac{dw}{w} = \frac{0.02 \text{ mm}}{0.6745(w)}.$$

$$\therefore \frac{dw_1}{w_1} = 0.145 \quad \frac{dw_2}{w_2} = 0.129 \quad \frac{dw_3}{w_3} = 0.114.$$



$$\begin{aligned}
 \text{Hence, } \frac{d\theta_0}{\theta_0} &= \frac{1}{\theta_0} * \frac{1}{\sqrt{2D}} \left(\frac{1}{2}\right) \frac{1}{(\omega_3^2 - 2\omega_2^2 + \omega_1^2)^{1/2}} * 2 \left( \frac{d\omega_3}{\omega_3} + \frac{2d\omega_2}{\omega_2} + \frac{d\omega_1}{\omega_1} \right) \\
 &= \frac{1}{\theta_0} * \frac{1}{2D} \left(\frac{1}{2}\right) \theta_0 * 2 (0.114 + 0.258 + 0.145) \\
 &= \frac{0.517}{200} \approx \underline{0.259\%}
 \end{aligned}$$

$$\text{and } \frac{d\omega_0}{\omega_0} = \left( \frac{d\lambda}{\lambda} + \frac{d\theta_0}{\theta_0} \right) \left( \frac{1}{\pi} \right) = \underline{0.082\%}$$

$$\begin{aligned}
 \text{and finally, } \frac{dz}{z} &= \frac{1}{z} d \left( \frac{1}{\theta_0} \cdot \sqrt{\omega_1^2 - \omega_0^2} \right) \\
 &= \frac{1}{z} \left( \frac{(-1)}{\theta_0^2} d\theta_0 \sqrt{\omega_1^2 - \omega_0^2} + \frac{1}{\theta_0} * \frac{1}{2\sqrt{\omega_1^2 - \omega_0^2}} * 2 \left( \frac{d\omega_1}{\omega_1} + \frac{d\omega_0}{\omega_0} \right) \right) \\
 &= \frac{1}{z} * \frac{1}{\theta_0} \left( \frac{d\theta_0}{\theta_0^2} (\omega_1^2 - \omega_0^2) + \frac{d\omega_1}{\omega_1} + \frac{d\omega_0}{\omega_0} \right) \\
 &= ((0.259)(0.0396) + 0.082 + 14.5) \% \\
 &= \underline{14.59\%}
 \end{aligned}$$

Conclusion: (i) In this way, various parameters of the He-Ne laser are measured.

(ii) Divergence of Laser ( $\theta_0$ ) =  $4.1 \times 10^{-3} \pm 0.259\%$

Beam Radius ( $\omega_0$ ) =  $49.13 \mu\text{m} \pm 0.082\%$

Distance of Beam Radii ( $z$ ) =  $48.537\text{mm} \pm 14.59\%$

Sources of Error: (i) Care should be taken to ensure that the laser beam, after reflection off of the detector does not enter into the laser apparatus or human eyes.  
(Precautions)

- (ii) Massive error is propagated due to  $\omega_1$ ,  $\omega_2$ , and  $\omega_3$ . Ensure that the Power v.s knife-edge distance graphs are sufficiently accurate.
- (iii) Ensure that the laser, knife-edge, and detector are perfectly aligned with respect to one another so that they form a straight line.